




PHLEBOLOGY


LEONARDO® DepthControl




Technology
meets Anatomy



The art of water absorption



Combined Efficacy – 1940 nm and 1470 nm joining together



Simplified action mechanism

Araujo et al Comparison of the effects of endovenous laser ablation at 1470 nm versus 1940 nm and different energy densities DOI: 10.1177/0268355518778488

1470 nm – Proven and effective

- Gold standard for safe & durable vein closure
- Also for large vein diameters (>15 mm)

1940 nm – Greater precision, gentler touch

- Lower penetration depth than 1470 nm – optimal control
- Particularly suitable for sensitive anatomical areas

Blending – It's the mix

Combining both wavelengths enables precise tailoring of the laser-tissue interaction to optimize recovery time, minimize side effects, and ensure reliable vein closure.

The energy delivery can be adapted to each vein's diameter, wall thickness, and proximity to nerves or other sensitive structures, allowing controlled thermal confinement. This dual-wavelength approach effectively limits lateral heat diffusion into surrounding tissue, thereby further reducing pain, hematomas, and nerve irritation.

Safety and comfort for your patients

- Less post-operative discomfort
- Faster recovery and less downtime
- Ultimate control





LEONARDO®

Model	LEONARDO® DUAL 20
REF	LL1470+1940nm20W
Wavelength	1470 nm ± 30 nm + 1940 nm ± 30 nm
Power	20 Watt (12 W 1470 nm + 8 W 1940 nm)
Fiber diameter	≥ 360 µm
Target beam	650 nm, red 4 mW, user controlled intensity
Treatment mode	CW, Pulse Mode, ELVeS® Signal, ELVeS® Segment
Pulse length / - break	0.01 – 60 sec / 0.01 – 60 sec
Power supply	100 – 240 VAC, 50/60 Hz, 450 VA
Dimensions (H × W × D)	approx. 28 cm × 37 cm × 9 cm
Weight	approx. 8.5 kg

All laser sets include three safety goggles, a foot switch, a door interlock, a power cable, and a manual in a carrying case.



Technical highlights

- Power: Up to 12 W (1470 nm) / 8 W (1940 nm)
- Modes: Single wavelength or combined use
- Display: Intuitive user interface with touchscreen
- Portable: Compact, lightweight design – ideal for clinics and practices

biolitec® worldwide

biolitec biomedical technology GmbH
Jena, Germany
Phone: +49 3641 519530

biolitec Schweiz GmbH
Wollerau, Switzerland
Phone: +41 55 5553020

biolitec España SL
Madrid, Spain
Phone: +34 91 9910857

biolitec Italia SRL
Milano, Italy
Phone: +39 02 84230633

biolitec Tıbbi Cihazları Ltd. Şti.
Istanbul, Turkey
Phone: +90 216 5747456

OOO biolitec Spb
Saint-Petersburg, Russia
Phone: +7 812 4493752

biolitec FZ LLC
Dubai, UAE
Phone: +971 44 298592

biolitec laser science and technology Shanghai Ltd.
Shanghai, China
Phone: +86 21 63088856

biolitec Sdn. Bhd.
Selangor, Malaysia
Phone: +60 3 55697158

biolitec India Private Ltd.
Bangalore, India
Phone: +91 9327511005

PT. Biolitec
Tangerang, Indonesia
Phone: +62 21 29557419

biolitec Korea Ltd.
Seoul, Republic of Korea
Phone: +82 2 7014707

Equipos Laser de Uso Médico y Fibra Óptica SA de CV
(Biolitec Mexico & Latin América)
Phone: +52 55 55731800

Biolitec Biotecnologia Comércio, Importação, Exportação LTDA
São Paulo, Brazil
Phone: +55 11 20938602



biolitec GmbH
Untere Viaduktgasse 6/9
A-1030 Wien
info@biolitec.com

Legal Manufacturer CeramOptec GmbH
Siemensstr. 44
53121 Bonn

Find your
local biolitec®
partner

