


LEONARDO[®] Duster Fibers

Leonardo Fibers



High-Quality, Durable & Stable Fibers for:

- Lithotripsy- Dusting,
Fragmentation & Popcorning
- Benign Prostatic Hyperplasia (BPH)
Enucleation & Vaporization
- Other Soft Tissue Treatment

LEONARDO® Duster ScopeSafe Fibers

REF	Product	Length [m]	Core $\varnothing$ [μ m]	OD $\varnothing$ [μ m]	Connector	NA
503200100	Leonardo® Duster ScopeSafe Fiber 150 μ m, IC	3.1	150	315	IC	0.26
503200105	Leonardo® Duster ScopeSafe Fiber 200 μ m, IC	3.1	200	400	IC	0.28
503200110	Leonardo® Duster ScopeSafe Fiber 272 μ m, IC	3.1	272	460	IC	0.22
503200115	Leonardo® Duster ScopeSafe Fiber 365 μ m, IC	3.1	365	550	IC	0.22
503200120	Leonardo® Duster ScopeSafe Fiber 550 μ m, IC	3.1	550	750	IC	0.22
503200125	Leonardo® Duster ScopeSafe Fiber 800 μ m, IC	3.1	800	1200	IC	0.22
503200130	Leonardo® Duster ScopeSafe Fiber 1000 μ m, IC	3.1	1000	1400	IC	0.22



**biolitec Holding
GmbH & Co KG**
Vienna, Austria
info@biolitec.de
www.biolitec.com

Legal Manufacturer
CeramOptec GmbH
Siemensstr. 44
53121 Bonn, Germany
Phone: +49 228 979670

Imprint
biolitec Holding GmbH & Co KG
Untere Viaduktgasse 6/9
A-1030 Wien
Phone: +43 1 3619 909 50
www.biolitec.com

All fibers are to be used only with the
1940 nm Leonardo® Duster Laser.