

Holmium und TFL Fasern



Besonders hochwertige und
exzellente Holmium und TFL Fasern




- Überragende Qualität
- Hohe Stabilität und Haltbarkeit
- Geeignet für diverse Lasersysteme



biolitec Holmium und TFL Fasern

biolitec fertigt diese qualitativ äußerst hochwertigen Laserfasern in der eigenen Produktion mit jahrzehntelangem Know-how. Ihre Klinik und Ihre Patienten profitieren von unseren High-Tech Fasern von herausragender Qualität.



Einmal Fasern



REF	Produktname	Länge [m]	Kern $\varnothing$ [μ m]/[Fr]	AD $\varnothing$ [μ m]/[Fr]	Stecker	NA
500200326	Megabeam HBF SF 230 - 253, OD 420 μ m	3.1	272	420	SMA905	0.28
500200327	Megabeam HBF SF 365 - 403, OD 550 μ m	3.1	365	550	SMA905	0.22/0.37
500200328	Megabeam HBF SF 550 - 603, OD 750 μ m	3.1	550	750	SMA905	0.22/0.37

Mehrfachnutzung



REF	Produktname	Länge [m]	Kern $\varnothing$ [μ m]/[Fr]	AD $\varnothing$ [μ m]/[Fr]	Stecker	NA
500200789	Megabeam RHBFSF 230 / 253 - 3	3.0	272	420	SMA905	0.22/0.36
500200790	Megabeam RHBFSF 365 / 400 - 3	3.0	365	550	SMA905	0.22/0.36
500200791	Megabeam RHBFSF 550 / 605 - 3	3.0	550	750	SMA905	0.22/0.36
500200792	Megabeam RHBFSF 945 / 1000 - 3	3.0	945	1400	SMA905	0.22/0.36

Side Fiber



REF	Produktname	Länge [m]	Kern $\varnothing$ [μ m]/[Fr]	AD $\varnothing$ [μ m]/[Fr]	Stecker	NA
500300451	Megabeam SF - 2100 - H	3.2	550	950	SMA905	0.22

Fasern sind kompatibel mit, aber nicht limitiert auf: StoneLight®, COHERENT®, Deka, LUMENIS®, Dornier MedTech und NewStarLaser Systeme.



**biolitec Holding
GmbH & Co KG**
Vienna, Austria
info@biolitec.de
www.biolitec.com

**Sales Department
biolitec biomedical
technology GmbH**
Jena, Germany
Phone: + 49 3641 519 53 0

Imprint
biolitec Holding GmbH & Co KG
Untere Viaduktgasse 6/9
A-1030 Wien
Phone: +43 1 3619 909 50
www.biolitec.com

All fibers are free of latex and DEHP.
Our fibers are single use products
(unless otherwise indicated) delivered
sterile for immediate use.

